

**Before you start**

See Medium & Reagents for JM8.F6 cells' protocol for all reagents

Thaw and pre-warm trypsin

Pre-warm PBS

Pre-warm required media

Pre-warm gelatin

Dish/Flask/Plate	Volume of Gelatin/PBS	Volume of Trypsin	Volume of Media
6 well plate	3ml	500 $\mu$ l	4.5ml
96 well plate	100 $\mu$ l	25 $\mu$ l	175 $\mu$ l
10cm Petri dish	8ml	1.5ml	8.5ml

**Things you'll need**

An aspirator & sterile filter-free tips/pipettes

Pipette aid & plastic pipettes

3 x liquid reservoirs

Multichannel pipette and 200 $\mu$ l Rainin filtered tips.

*We have tried several brands of tips, and find that the bore of the Rainin 200 $\mu$ l tip is optimal for dissociating cells.*

**For passaging cells in 96-well plates (Falcon #353916) or 10cm<sup>2</sup> Petri dishes (Corning #430167):**

*Example: Passing a confluent 10cm Petri dish into 4 daughter Petri dishes.*

1. Remove the parent 10cm Petri dish from the incubator and observe under the microscope. Check cells are healthy and confluent
2. Aspirate old media
3. Wash with 10 ml pre-warmed PBS  
*Add PBS to the side of the dish, and slowly tilt dish to gently wash the cells.*
4. Add 1.5 ml pre-warmed trypsin.  
*Gently swirl the dish to cover all cells with trypsin*
5. Incubate at 37° for 15 minutes
6. While cells are incubating, remove medium from fresh feeder plates and add fresh, pre-warmed medium
7. After incubation, gently swirl the plate again and add 8.5 ml medium to inactivate the trypsin.  
*Pipette up and down gently 3 to 4 times to disperse cells*
8. Transfer 2.5 ml of the cell suspension to each of the four fresh feeder plates.  
*Swirl the plate to distribute the cells evenly across the plate*
9. Incubate the plates in an incubator at 37°C with 5% CO<sub>2</sub>.
10. Check cells the following day and media change to remove all traces of trypsin and dead cells.

*Example: Passing a confluent 96 well plate into 4 daughter plates.*

1. Coat four fresh, sterile 96-well tissue culture plates with 100  $\mu$ l gelatin.  
*Incubate at room temperature for >10 minutes.*
2. Remove the parent 96-well plate from the incubator and observe under the microscope. Check cells are healthy and confluent
3. Aspirate old media
4. Wash each well with 100  $\mu$ l pre-warmed PBS  
*Add PBS to the side of the wells*
5. Add 25  $\mu$ l pre-warmed trypsin.  
*Gently tap the plate to ensure all cells are covered with trypsin*

6. Incubate at 37° incubator for 15 minutes
7. While cells are incubating, remove gelatin from the daughter 96-well plates and add 150 µl fresh, pre-warmed medium to each well
8. After incubation, gently tap the plate again and add 175 µl medium to inactivate the trypsin.  
*Gently pipette up and down 3 to 4 times to disperse cells.*
9. Gently transfer 50 µl of the cell suspension to each of the four daughter 96-well plates.
10. Incubate the plates in an incubator at 37°C with 5% CO<sub>2</sub>.
11. Check cells the following day and media change to remove all traces of trypsin and dead cells

---

**12. For splitting 96-well plates:**

*Example: Splitting a 96-well plate 1:4.*

- 1) Gelatinise 4 X 96-well plates with 100µl 0.1% gelatin.
- 2) Remove the 96-well plate from the incubator and observe under the microscope. Check cells are healthy and confluent.
- 3) Using multichannel aspirator, remove old media from wells.
- 4) Using electronic multichannel pipette, add 100µl of PBS to each well.
- 5) Aspirate off PBS.
- 6) Add 25µl 2 x Trypsin + glucose.
- 7) Gently tap plate and incubate for 8-10 minutes.
- 8) While cells are incubating, remove gelatin using aspirator and add 150µl of pre-warmed media to each well.
- 9) After incubation, gently tap the plate again, and check that cells have lifted off the plate.
- 10) Add 175µl of media to the cells.
- 11) Use the Eppendorf multichannel pipette and the 200µl Rainin tips for the split:
  - i. Use a fresh tip for each well i.e. a whole box of 96 tips will be required
  - ii. Gently aspirate the first column of the master plate, up and down about 5 times, then transfer 50µl of cell suspension to the first column of each of the new plates, gently mixing the cell suspension in the new media.
  - iii. Discard tips, and repeat for the next 11 columns.
- 12) Label plates and incubate.
- 13) Check cells the following day and media change to remove all traces of trypsin and dead cells.