

Before you start

See 'Medium & Reagents for JM8.N4 cells' protocol for all reagents

Thaw and pre-warm trypsin

Pre-warm PBS

Pre-warm required media

Pre-warm gelatin

| Dish/Flask/Plate | Volume of Gelatin/PBS | Volume of Trypsin | Volume of Media |
|------------------|-----------------------|-------------------|-----------------|
| 6 well plate | 3ml | 500 μ l | 4.5ml |
| 96 well plate | 100 μ l | 25 μ l | 175 μ l |
| 10cm Petri dish | 8ml | 1.5ml | 8.5ml |

Things you'll need

An aspirator & sterile filter-free tips/pipettes

Pipette aid & plastic pipettes

3 x liquid reservoirs

Multichannel pipette and 200 μ l Rainin filtered tips.

We have tried several brands of tips, and find that the bore of the Rainin 200 μ l tip is optimal for dissociating cells.

For passaging cells in 96-well plates (Falcon #353916) or flasks (Corning #430167):

JM8.N4 cells should be passaged every 48 hours at a 1:6 dilution

1. Remove the parent *flask* from the incubator and observe under the microscope. Check cells are healthy and confluent
2. Aspirate old media
3. Wash with pre-warmed PBS
Add PBS to the side of the flask, and slowly tilt to gently wash the cells.
4. Add pre-warmed trypsin.
Gently swirl the flask to cover all cells with trypsin
5. Incubate at 37° for 10 minutes
6. While cells are incubating, coat fresh flasks with 0.1% gelatine (1 minute), aspirate and add pre-warmed media
7. After incubation, gently swirl the flask again and add 5 volume of medium to inactivate the trypsin.
Pipette up and down gently 3 to 4 times to disperse cells
8. Transfer 1/6th of cell suspension to each of the fresh flasks.
9. Incubate the plates in an incubator at 37°C with 5% CO₂.
10. Check cells the following day and media change to remove all traces of trypsin and dead cells.

Example: Passing a confluent 96 well plate into 4 daughter plates.

1. Coat four fresh, sterile 96-well tissue culture plates with 100 μ l gelatin.
Incubate at room temperature for >10 minutes.
2. Remove the parent 96-well plate from the incubator and observe under the microscope. Check cells are healthy and confluent
3. Aspirate old media
4. Wash each well with 100 μ l pre-warmed PBS
Add PBS to the side of the wells
5. Add 25 μ l pre-warmed trypsin.
Gently tap the plate to ensure all cells are covered with trypsin
6. Incubate at 37° incubator for 10 minutes

7. While cells are incubating, remove gelatin from the daughter 96-well plates and add 150 μ l fresh, pre-warmed medium to each well
8. After incubation, gently tap the plate again and add 175 μ l medium to inactivate the trypsin.
Gently pipette up and down 3 to 4 times to disperse cells.
9. Gently transfer 50 μ l of the cell suspension to each of the four daughter 96-well plates.
10. Incubate the plates in an incubator at 37°C with 5% CO₂.
11. Check cells the following day and media change to remove all traces of trypsin and dead cells

12. For splitting 96-well plates:

Example: Splitting a 96-well plate 1:4.

- 1) Gelatinise 4 X 96-well plates with 100µl 0.1% gelatin.
- 2) Remove the 96-well plate from the incubator and observe under the microscope. Check cells are healthy and confluent.
- 3) Using multichannel aspirator, remove old media from wells.
- 4) Using electronic multichannel pipette, add 100µl of PBS to each well.
- 5) Aspirate off PBS.
- 6) Add 25µl 2 x Trypsin + glucose.
- 7) Gently tap plate and incubate for 8-10 minutes.
- 8) While cells are incubating, remove gelatin using aspirator and add 150µl of pre-warmed media to each well.
- 9) After incubation, gently tap the plate again, and check that cells have lifted off the plate.
- 10) Add 175µl of media to the cells.
- 11) Use the Eppendorf multichannel pipette and the 200µl Rainin tips for the split:
 - i. Use a fresh tip for each well i.e. a whole box of 96 tips will be required
 - ii. Gently aspirate the first column of the master plate, up and down about 5 times, then transfer 50µl of cell suspension to the first column of each of the new plates, gently mixing the cell suspension in the new media.
 - iii. Discard tips, and repeat for the next 11 columns.
- 12) Label plates and incubate.
- 13) Check cells the following day and media change to remove all traces of trypsin and dead cells.